

Power Mobility

Complex Rehab Vs. Basic



- Only 7% of Medicare Power WCs
- Intended for Progressive Diagnoses
- Advanced Electronics and Controls
- Provides Positioning
- Accommodates Deformity
- Provides Pressure Management
- Offers Ventilator Accommodation



- 93% of Medicare Power WCS
- Intended for Ambulatory Limitations
- Basic Joystick Drive ONLY
- Provides NO Positioning
- NO Deformity Accommodation
- Provides NO Pressure Management
- NO Ventilator Accommodation

Manual Wheelchairs

Complex Rehab Vs. Basic



- Only 6% of Medicare Manual WCs
- Intended for long-term use
- High Adjustability
- Provides Positioning
- Accommodates Deformity
- Provides Pressure Management



- 94% of Medicare Manual WCS
- Intended for short-term use
- Minimal to Zero Adjustability
- Provides NO Positioning
- NO Deformity Accommodation
- Provides NO Pressure Management